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Milli-Q® Direct 8/16 系統

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
步驟	動作	圖
4	取水至少 10 分鐘。	 <p>A screenshot of a control panel display. At the top, it says 'READY'. Below that is the date and time: '02 jul. 2009 13:55'. There are three menu options: 'Menu →', 'Tank: Standby →', and '80.0 % Volume →'. Below these are several water quality parameters: 'Perm C: 6.0 µS/cm TC', 'MQ Res: 18.2 MΩcm TC', and 'TOC: 4 ppb'.</p>
5	將 POD 開關向下按，然後鬆開，以停止取水。 使系統停留在 READY 模式。	 <p>A screenshot of a control panel display, identical to the one in step 4. It shows 'READY' status, the date and time '02 jul. 2009 13:55', menu options 'Menu →', 'Tank: Standby →', and '80.0 % Volume →', and water quality parameters: 'Perm C: 6.0 µS/cm TC', 'MQ Res: 18.2 MΩcm TC', and 'TOC: 4 ppb'.</p>

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
圖 1	圖 2																							
<table border="1"> <tr><td>STANDBY MENU</td></tr> <tr><td>Maintenance →</td></tr> <tr><td>Sanitise/Clean →</td></tr> <tr><td>Suitability Tests →</td></tr> <tr><td>Language →</td></tr> <tr><td>Manager Menu →</td></tr> </table>	STANDBY MENU	Maintenance →	Sanitise/Clean →	Suitability Tests →	Language →	Manager Menu →	<table border="1"> <tr><td>MAINTENANCE</td></tr> <tr><td>Install Pretreatment →</td></tr> <tr><td>Clean Strainer →</td></tr> <tr><td>Install Progard →</td></tr> <tr><td>Install new RO →</td></tr> <tr><td>Install UV 254 nm Lamp →</td></tr> <tr><td>Install UV 185 nm Lamp →</td></tr> <tr><td>Install Q-Pak →</td></tr> </table>	MAINTENANCE	Install Pretreatment →	Clean Strainer →	Install Progard →	Install new RO →	Install UV 254 nm Lamp →	Install UV 185 nm Lamp →	Install Q-Pak →	<table border="1"> <tr><td>MAINTENANCE</td></tr> <tr><td>Install new RO →</td></tr> <tr><td>Install UV 254 nm Lamp →</td></tr> <tr><td>Install UV 185 nm Lamp →</td></tr> <tr><td>Install Q-Pak →</td></tr> <tr><td>Install POD Pak →</td></tr> <tr><td>Install ASM UV lamp →</td></tr> <tr><td>Reset RO CL2 CLEANING →</td></tr> </table>	MAINTENANCE	Install new RO →	Install UV 254 nm Lamp →	Install UV 185 nm Lamp →	Install Q-Pak →	Install POD Pak →	Install ASM UV lamp →	Reset RO CL2 CLEANING →
STANDBY MENU																								
Maintenance →																								
Sanitise/Clean →																								
Suitability Tests →																								
Language →																								
Manager Menu →																								
MAINTENANCE																								
Install Pretreatment →																								
Clean Strainer →																								
Install Progard →																								
Install new RO →																								
Install UV 254 nm Lamp →																								
Install UV 185 nm Lamp →																								
Install Q-Pak →																								
MAINTENANCE																								
Install new RO →																								
Install UV 254 nm Lamp →																								
Install UV 185 nm Lamp →																								
Install Q-Pak →																								
Install POD Pak →																								
Install ASM UV lamp →																								
Reset RO CL2 CLEANING →																								

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
圖 1	圖 2
 <p>STANDBY MENU Maintenance → Sanitise/Clean → Suitability Tests → Language → Manager Menu →</p>	 <p>SANITISE / CLEAN RO CL2 Cleaning → RO pH Cleaning → RO Cleaning → System Cleaning →</p>

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
圖 1	圖 2
 <p>STANDBY MENU Maintenance → Sanitise/Clean → Suitability Tests → Language → Manager Menu →</p>	 <p>SUITABILITY TESTS Res Suitability Test → Temp Suitability Test →</p>

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圖 1	圖 2
 <p>STANDBY MENU Maintenance → Sanitise/Clean → Suitability Tests → Language → Manager Menu →</p>	 <p>LANGUAGE Chinese English ✓ French German Italian Japanese Portuguese</p>

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
圖 1	圖 2																									
<table border="1"> <tr><td>MANAGER MENU</td></tr> <tr><td>Change ID and Password →</td></tr> <tr><td>Date and Time →</td></tr> <tr><td>Set Points →</td></tr> <tr><td>Units →</td></tr> <tr><td>Setup →</td></tr> <tr><td>User Parameters →</td></tr> <tr><td>History →</td></tr> </table>	MANAGER MENU	Change ID and Password →	Date and Time →	Set Points →	Units →	Setup →	User Parameters →	History →	<table border="1"> <tr><td>SET POINTS</td></tr> <tr><td>Pretreatment →</td></tr> <tr><td>Strainer Frequency →</td></tr> <tr><td>Tap Feed Cond →</td></tr> <tr><td>RO Rejection →</td></tr> <tr><td>RO CL2 Cleaning →</td></tr> <tr><td>Permeate Cond →</td></tr> <tr><td>Tank Refill →</td></tr> </table>	SET POINTS	Pretreatment →	Strainer Frequency →	Tap Feed Cond →	RO Rejection →	RO CL2 Cleaning →	Permeate Cond →	Tank Refill →	<table border="1"> <tr><td>SET POINTS</td></tr> <tr><td>Milli-Q Product Res →</td></tr> <tr><td>Milli-Q Product TOC →</td></tr> <tr><td>Millipak →</td></tr> <tr><td>BioPak →</td></tr> <tr><td>EDS-Pak →</td></tr> <tr><td>Pod Pak A →</td></tr> <tr><td>Pod Pak B →</td></tr> </table>	SET POINTS	Milli-Q Product Res →	Milli-Q Product TOC →	Millipak →	BioPak →	EDS-Pak →	Pod Pak A →	Pod Pak B →
MANAGER MENU																										
Change ID and Password →																										
Date and Time →																										
Set Points →																										
Units →																										
Setup →																										
User Parameters →																										
History →																										
SET POINTS																										
Pretreatment →																										
Strainer Frequency →																										
Tap Feed Cond →																										
RO Rejection →																										
RO CL2 Cleaning →																										
Permeate Cond →																										
Tank Refill →																										
SET POINTS																										
Milli-Q Product Res →																										
Milli-Q Product TOC →																										
Millipak →																										
BioPak →																										
EDS-Pak →																										
Pod Pak A →																										
Pod Pak B →																										

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圖 1	圖 2
 <p>MANAGER MENU Change ID and Password → Date and Time → Set Points → Units → Setup → User Parameters → History →</p>	 <p>UNITS Pressure → Milli-Q Product → Tank Volume →</p>

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圖 1	圖 2																								
<table border="1"> <tr><td>MANAGER MENU</td></tr> <tr><td>Change ID and Password →</td></tr> <tr><td>Date and Time →</td></tr> <tr><td>Set Points →</td></tr> <tr><td>Units →</td></tr> <tr><td>Setup →</td></tr> <tr><td>User Parameters →</td></tr> <tr><td>History →</td></tr> </table>	MANAGER MENU	Change ID and Password →	Date and Time →	Set Points →	Units →	Setup →	User Parameters →	History →	<table border="1"> <tr><td>SETUP</td></tr> <tr><td>Install Date →</td></tr> <tr><td>Buzzer →</td></tr> <tr><td>MQ Recirc Mode →</td></tr> <tr><td>POD Flow Stop →</td></tr> <tr><td>Temp Comp Mode →</td></tr> <tr><td>Flow Calibration →</td></tr> <tr><td>UV 254 nm Activation →</td></tr> </table>	SETUP	Install Date →	Buzzer →	MQ Recirc Mode →	POD Flow Stop →	Temp Comp Mode →	Flow Calibration →	UV 254 nm Activation →	<table border="1"> <tr><td>SETUP</td></tr> <tr><td>Flow Calibration →</td></tr> <tr><td>UV 254 nm Activation →</td></tr> <tr><td>UV 185 nm Activation →</td></tr> <tr><td>ASM UV Lamp Schedule →</td></tr> <tr><td>Network Settings →</td></tr> <tr><td>TOC Curve Check Time →</td></tr> </table>	SETUP	Flow Calibration →	UV 254 nm Activation →	UV 185 nm Activation →	ASM UV Lamp Schedule →	Network Settings →	TOC Curve Check Time →
MANAGER MENU																									
Change ID and Password →																									
Date and Time →																									
Set Points →																									
Units →																									
Setup →																									
User Parameters →																									
History →																									
SETUP																									
Install Date →																									
Buzzer →																									
MQ Recirc Mode →																									
POD Flow Stop →																									
Temp Comp Mode →																									
Flow Calibration →																									
UV 254 nm Activation →																									
SETUP																									
Flow Calibration →																									
UV 254 nm Activation →																									
UV 185 nm Activation →																									
ASM UV Lamp Schedule →																									
Network Settings →																									
TOC Curve Check Time →																									

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圖 1	圖 2																
<table border="1"> <tr><td>MANAGER MENU</td></tr> <tr><td>Change ID and Password →</td></tr> <tr><td>Date and Time →</td></tr> <tr><td>Set Points →</td></tr> <tr><td>Units →</td></tr> <tr><td>Setup →</td></tr> <tr><td>User Parameters →</td></tr> <tr><td>History →</td></tr> </table>	MANAGER MENU	Change ID and Password →	Date and Time →	Set Points →	Units →	Setup →	User Parameters →	History →	<table border="1"> <tr><td>USER PARAMETERS</td></tr> <tr><td>Company Name →</td></tr> <tr><td>Department Name →</td></tr> <tr><td>Address →</td></tr> <tr><td>Postal Code →</td></tr> <tr><td>City →</td></tr> <tr><td>Country →</td></tr> <tr><td>Email →</td></tr> </table>	USER PARAMETERS	Company Name →	Department Name →	Address →	Postal Code →	City →	Country →	Email →
MANAGER MENU																	
Change ID and Password →																	
Date and Time →																	
Set Points →																	
Units →																	
Setup →																	
User Parameters →																	
History →																	
USER PARAMETERS																	
Company Name →																	
Department Name →																	
Address →																	
Postal Code →																	
City →																	
Country →																	
Email →																	

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圖 1	圖 2														
<table border="1"> <tr><td>MANAGER MENU</td></tr> <tr><td>Date and Time →</td></tr> <tr><td>Set Points →</td></tr> <tr><td>Units →</td></tr> <tr><td>Setup →</td></tr> <tr><td>User Parameters →</td></tr> <tr><td>History →</td></tr> <tr><td>Lab closed →</td></tr> </table>	MANAGER MENU	Date and Time →	Set Points →	Units →	Setup →	User Parameters →	History →	Lab closed →	<table border="1"> <tr><td>HISTORY</td></tr> <tr><td>History Summary →</td></tr> <tr><td>Print System History →</td></tr> <tr><td>Print RO History →</td></tr> <tr><td>Print Milli-Q History →</td></tr> <tr><td>Print Options →</td></tr> </table>	HISTORY	History Summary →	Print System History →	Print RO History →	Print Milli-Q History →	Print Options →
MANAGER MENU															
Date and Time →															
Set Points →															
Units →															
Setup →															
User Parameters →															
History →															
Lab closed →															
HISTORY															
History Summary →															
Print System History →															
Print RO History →															
Print Milli-Q History →															
Print Options →															

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


圖 1	圖 2
 <p>READY MENU Water Quality → Print Menu → View Operation → Consumables Status → Call Millipore → Service Tracking → Information →</p>	 <p>WATER QUALITY Permeate Water Quality → Tank Level: 88.0 % Milli-Q Water Quality →</p>

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

圖 1	圖 2
 <p>READY MENU Water Quality → Print Menu → View Operation → Consumables Status → Call Millipore → Service Tracking → Information →</p>	 <p>PRINT MENU MQ Instant Quality →</p>

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圖 1	圖 2
 <p>READY MENU Water Quality → Print Menu → View Operation → Consumables Status → Call Millipore → Service Tracking → Information →</p>	 <p>VIEW OPERATION System Operation → System Alerts → System Alarms → System Measures →</p>

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圖 1	圖 2																
<table border="1"> <tr><td>READY MENU</td></tr> <tr><td>Water Quality →</td></tr> <tr><td>Print Menu →</td></tr> <tr><td>View Operation →</td></tr> <tr><td>Consumables Status →</td></tr> <tr><td>Call Millipore →</td></tr> <tr><td>Service Tracking →</td></tr> <tr><td>Information →</td></tr> </table>	READY MENU	Water Quality →	Print Menu →	View Operation →	Consumables Status →	Call Millipore →	Service Tracking →	Information →	<table border="1"> <tr><td>CONSUMABLES STATUS</td></tr> <tr><td>Pretreatment →</td></tr> <tr><td>Progard →</td></tr> <tr><td>UV 254 nm Lamp →</td></tr> <tr><td>ASM UV Lamp →</td></tr> <tr><td>UV 185 nm Lamp →</td></tr> <tr><td>Q-Pak →</td></tr> <tr><td>POD Pak →</td></tr> </table>	CONSUMABLES STATUS	Pretreatment →	Progard →	UV 254 nm Lamp →	ASM UV Lamp →	UV 185 nm Lamp →	Q-Pak →	POD Pak →
READY MENU																	
Water Quality →																	
Print Menu →																	
View Operation →																	
Consumables Status →																	
Call Millipore →																	
Service Tracking →																	
Information →																	
CONSUMABLES STATUS																	
Pretreatment →																	
Progard →																	
UV 254 nm Lamp →																	
ASM UV Lamp →																	
UV 185 nm Lamp →																	
Q-Pak →																	
POD Pak →																	

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
圖 1	圖 2													
<table border="1"> <tr><td>READY MENU</td></tr> <tr><td>Water Quality →</td></tr> <tr><td>Print Menu →</td></tr> <tr><td>View Operation →</td></tr> <tr><td>Consumables Status →</td></tr> <tr><td>Call Millipore →</td></tr> <tr><td>Service Tracking →</td></tr> <tr><td>Information →</td></tr> </table>	READY MENU	Water Quality →	Print Menu →	View Operation →	Consumables Status →	Call Millipore →	Service Tracking →	Information →	<table border="1"> <tr><td>CALL MILLIPORE</td></tr> <tr><td>Application Specialist →</td></tr> <tr><td>Service Engineer →</td></tr> <tr><td>Tech Service →</td></tr> <tr><td>Other →</td></tr> </table>	CALL MILLIPORE	Application Specialist →	Service Engineer →	Tech Service →	Other →
READY MENU														
Water Quality →														
Print Menu →														
View Operation →														
Consumables Status →														
Call Millipore →														
Service Tracking →														
Information →														
CALL MILLIPORE														
Application Specialist →														
Service Engineer →														
Tech Service →														
Other →														

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
圖 1	圖 2
 <p>READY MENU Water Quality → Print Menu → View Operation → Consumables Status → Call Millipore → Service Tracking → Information →</p>	 <p>SERVICE TRACKING Installation → Repair → Service Contract → Contract Expires → Next Service → Next Calibration → Next Qualification →</p>

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
圖 1	圖 2
 <p>READY MENU Water Quality → Print Menu → View Operation → Consumables Status → Call Millipore → Service Tracking → Information →</p>	 <p>INFORMATION Flow Schematic → Version → System Information →</p>

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圖 1	圖 2
 <p>READY MENU Water Quality → View Operation → Consumables Status → Call Millipore → Service Tracking → InFormation → TOC Curve Check →</p>	 <p>TOC CURVE CHECK Press ✓ to start TOC curve check operation. Press ← to exit.</p>

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© 2009 密理博公司。法國印刷。版權所有。未經出版商書面許可，不得以任何方式複製本書的全部或部分內容。